Replenish™

Multi-Purpose Disinfecting Solution
Διάλυμα Απολύμανσης Πολλαπλών Χρήσεων
Solução Multiusos Desinfetante
Solución Desinfectante Multi-acción

INSTRUCTIONS FOR USE

EN This Instructions For Use leaflet contains important product use and safety information. Please read carefully and keep this for future use.

PRODUCT NAME:

OPTI-FREE™ REPLENISH™ Multi-Purpose Disinfecting Solution and Alcon Contact Lens Case (Flat)

DESCRIPTION:

OPTI-FREE™ REPLENISH™ Multi-Purpose Disinfecting Solution and Alcon Contact Lens Case (Flat). OPTI-FREE™ REPLENISH™ Solution is a sterile, buffered, aqueous solution containing sodium citrate, sodium chloride, sodium borate, propylene glycol, TearGlyde™ (proprietary dual action reconditioning system (TETRONIC™ 1304*, nonanol ethylenediaminetetraacetic acid) with POLYQUAD™ (polyquaternium-1) 0.001% and ALDOX™ (myristamidopropyl dimethylamine) 0.0005% preservatives, and purified water. The Alcon Contact Lens Case (Flat) is an accessory to OPTI-FREE™ REPLENISH™ Solution. The lens case is a flat, side-by-side well design and consists of two closures and a base with two lens well stations. OPTI-FREE™ REPLENISH™ Solution and Alcon Contact Lens Case should always be used together

CLINICAL BENEFITS:

When used as directed, OPTI-FREE™ REPLENISH™ Solution:

- **Cleans** by loosening and removing accumulations of protein and other deposits and debris from your lenses during disinfection and storage.
- **Disinfects**: Biocompatible, dual disinfectants POLYQUAD™ and ALDOX™ have activity against harmful microorganisms.
- **Rinses** and stores soft contact lenses (including silicone hydrogel lenses).
- **Reconditions** the surface of the lens thereby improving lens surface, lowering lens deposits and improved lens wettability.
- **Stores** lenses in an unopened lens case for up to 30 days after disinfection.
- Provides contact lens comfort.

INTENDED PURPOSE/USE:

OPTI-FREE™ REPLENISH™ Solution is intended for cleaning, disinfection, and storage of soft contact lenses (including silicone hydrogel lenses).

INTENDED USES/TARGET POPULATION:

Contact lens wearers (7 years and older) using soft contact lenses (including silicone hydrogel lenses), who are capable of understanding and following the directions and instructions for use (or have a caregiver who can act on their behalf), and eye care professionals.

PERFORMANCE CHARACTERISTICS:

Based on the total body of evidence from clinical investigations, published clinical literature and clinical experience (postmarket surveillance), OPTI-FREE™ REPLENISH™ Solution performs as intended for cleaning, disinfection, and storage of soft contact lenses (including silicone hydrogel lenses).

DIRECTIONS FOR USE:

General Directions:

- Remove imprinted seal around bottleneck completely before use.
- Always wash and rinse your hands thoroughly and dry with a clean, lint free towel before handling your lenses.
- Clean, rinse and disinfect your lenses each time you remove them.
- Always handle the same lens, the right or the left, first in order to avoid mix-ups.
- When using a new Alcon Contact Lens Case, rinse the case with OPTI-FREE™ REPLENISH™ Solution before use to remove any debris.

Follow These Directions for Use:

- **Thoroughly wet** each side of the lens with OPTI-FREE™ REPLENISH™ Solution.
- **Rub** each side of the lens for 10 seconds (for a total rubbing time of 20 seconds).
- **Rinse** each side of the lens for 10 seconds with a steady stream of OPTI-FREE™ REPLENISH™ Solution (for a total rinsing time of 20 seconds).
- **Fill** your Alcon Contact Lens Case with fresh OPTI-FREE™ REPLENISH™ Solution. **Store** lenses in the closed lens case overnight or at least 6 hours. After soaking, lenses are disinfected and ready to wear.
- If any debris remains on contact lenses, rinse with OPTI-FREE™ REPLENISH™ Solution prior to lens insertion. You may leave your lenses in the unopened lens case containing OPTI-FREE™ REPLENISH™ Solution for up to 30 days. After this time, your lenses must be cleaned and disinfected with fresh OPTI-FREE™ REPLENISH™ Solution prior to wear.

Cleaning of lens case:

- After use, always empty and clean the Alcon Contact Lens Case by rinsing with fresh OPTI-FREE™ REPLENISH™ Solution. Do not wash the lens case with water.
- Before allowing the lens case to air dry, be sure that no residual solution remains in the case.

CONTRAINDICATIONS (Reasons Not to Use):

Do not use this product if you are allergic to any of the ingredients.

WARNINGS:

It is essential that you follow your eye care professional's directions and all labeling instructions for proper use and care of your lenses and lens care products, including the lens case. Do not reuse or "top off" solution in the lens case since this may reduce effective lens disinfection. Eye problems, including an infection (serious microbial infection) or corneal ulcers, can develop rapidly and lead to loss of vision. If you experience persistent eye discomfort, eye pain, burning and stinging, dry eye sensation, excessive tearing, redness of the eye (irritation, allergic response), or vision changes (decreased vision), IMMEDIATELY remove your lenses and consult your eye care professional as the problem could become more serious.

PRECAUTIONS:

- Always use fresh OPTI-FREE™ REPLENISH™ Solution. Never reuse the solution in your lens storage case.
- Always use the new Alcon Contact Lens Case (Flat) provided with each purchase of OPTI-FREE™ REPLENISH™ Solution.
- To help prevent eye infections, replace your lens case frequently.
- TAMPER EVIDENT: Do not use if the security seal is damaged or missing prior to first opening.
- Do not use if package or bottle (sterile packaging) or lens case is damaged.
- Never use water, saline solution, or rewetting drops to disinfect your lenses.
- Do not swallow. If swallowed, get medical help as abdominal discomfort (gastrointestinal discomfort) may occur.
- To avoid contamination, never touch dropper tip of the container to any surface.
- Close the flip-top cap after using.
- Store between 15°C to 30°C.
- Keep the container tightly closed when not in use.
- This product contains a level of sodium borate that is safe when used according to label instructions. Sodium borate has a potential to cause reproductive toxicity.
- Use before the expiration date marked on the product.
- Once the bottle is opened, discard any remaining solution after six months.
- Adult supervision required. Keep out of the reach of children.
- Not for use with heat (thermal) disinfection.
- Do not change lens care solutions without consulting your eye care professional.

RESIDUAL RISKS AND UNDESIRABLE SIDE EFFECTS:

If you experience persistent eye discomfort, eye pain, burning and stinging, dry eye sensation, excessive tearing, gastrointestinal (abdominal) discomfort, redness of the eye (irritation, allergic response/hypersensitivity), or vision changes (blurred vision, decreased vision), immediately remove your lenses. Consult your eye care professional to identify the problem and seek treatment to avoid serious injury. An eye infection (serious microbial infection) may be present and if untreated, could lead to loss of vision (permanent decreased vision).

These signs and symptoms may be associated with multiple conditions, including Inflammation and Toxicity response (general and systemic).

Serious Incident Reporting:

Any serious incident related to the use of this medical device should be reported to Alcon Laboratories, Inc.:

Contact the local country office or your Alcon distributor.

Email: ga.complaints@alcon.com

Website: <https://www.alcon.com/contact-us>

These serious incidents should also be reported to the competent authority for medical devices of your country.

HOW SUPPLIED:

OPTI-FREE™ REPLENISH™ Solution, which is sterilized by filtration, is supplied in plastic bottles containing fill sizes ranging from 60 mL to 420 mL.

Alcon Lens Case included.

STERILE

A Sterilized using aseptic processing techniques

GLOSSARY TABLE FOR SYMBOLS

Symbol	Definitions
	Caution
	Consult instructions for use or consult electronic instructions for use
	Do not use if package is damaged and consult instructions for use
	Date of manufacture
	Contains hazardous substances
	Medical device
	Unique device identifier
	Single sterile barrier system
	Packaging waste license sign

Manufacturer
Alcon Laboratories, Inc.
6201 South Freeway
Fort Worth, TX 76134-2099, USA

© 2026 Alcon Inc.

TETRONIC is a trademark of BASF.

www.ifu.alcon.com

Date of issue: 2026-03

EC REP Authorized Representative in the European Union
Alcon Laboratories Belgium
Lichterveld 3
2870 Puurs-Sint-Amands
Belgium

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