02b Ophthalmic antipruritic solution indicated for the relief of symptoms of mild or moderate eye dryness. **PRESERVATIVE FREE. DESCRIPTION:** **LacriShield Oxy-Dry 0.2g** is a lubricating, hypotonic ophthalmic preservative-free ophthalmic solution with a combination of polyvinylpyrrolidone (sodium salt) and polyethylene glycol essential electrolytes, intended to relieve the symptoms of mild or moderate ocular dryness and compatible with all contact lens materials. **USAGE:** **LacriShield Oxy-Dry 0.2%** is indicated for:

- Extended relief of the symptoms of mild or moderate ocular dryness and other discomforts, such as stinging, itching, burning, irritation, or foreign body sensations associated with dry eyes.
- Relief of dryness, irritation, and burning associated with cigarette smoking, air conditioning, and prolonged use of computers and contact lens wear, among others.
- Lubricating solution for contact lenses.

INSTRUCTIONS FOR USE: 1. Wash, rinse and dry hands before using this product. 2. Take a bottle of **LacriShield Oxy-Dry 0.2%**, pull the safety seal antilockwise, and remove the protective cap. 3. Direct use: Apply one or two drops only on the eye and blink several times. 4. Indirect use: Apply one or two drops only on the external surface of each lens during wear or before taking them out, or apply one or two drops into the interior part of each lens before placing on the eye. 5. After each administration, place the cap back on. 6. Use the product, according to your needs, up to 8 times a day. Last revision date: 01/2017. **CAUTION:** Do not use the product if you are allergic to any of its ingredients. • Never use tap water to rinse the lens case or the lenses. • If you experience ocular irritation or an adverse reaction, please consult a specialist. • Once the bottle is opened, use the product within 90 days. • Ensure that the tip of the bottle does not come into contact with the eyes, or any other surface, to avoid contamination. • Do not use the product if the bottle is damaged. • Store the product at room temperature. • Do not use the product after the expiry date. • Keep the bottle tightly closed after each use. • Do not use the product in combination with a pharmacological treatment or any other product that is applied on the ocular surface, without consulting your doctor/eye care professional. **INGREDIENTS:** Polyvinylpyrrolidone (sodium salt), Polyethylene Glycol 0.20 g, Glycerine, Boric Acid, Potassium Chloride, Magnesium Chloride, Potassium Chloride, Calcium Chloride, Povidone, Borax Decahydrate and Purified Water 100 ml.

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