AVIZOR

Unica[®]
sensitive

WITH
hyaluronic acid

for sensitive eyes

e

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gb

ALL IN ONE SOLUTION FOR USERS WITH SENSITIVE EYES. FOR ALL TYPES OF SOFT CONTACT LENSES (including the new generation of silicone hydrogel-lenses).

UNICA SENSITIVE has been especially formulated to provide extra comfort to even the most sensitive contact lens users.

Functions: Cleans and eliminates proteins, disinfects, stores, rinses, wets and lubricate.

INSTRUCTIONS FOR USE:

1. Before handling your lenses, wash your hands with soap and water.
 2. Place your lens on the palm of your hand and rub it carefully with a few drops of UNICA SENSITIVE SOLUTION, on both sides for a few seconds.
 3. Rinse your lenses with UNICA SENSITIVE SOLUTION.
 4. Place your lenses in the lens case, and be sure to cover these with UNICA SENSITIVE SOLUTION.
 5. Keep lenses at least 4 hours in UNICA SENSITIVE SOLUTION (or overnight) to grant an adequate disinfection.
 6. Once the disinfection process is done, you may insert your lenses directly to your eye or you may rinse them previously with UNICA SENSITIVE SOLUTION.
 7. Keep your lens case clean and be sure to replace the UNICA SENSITIVE SOLUTION in your lens case daily. Do not use tap water to clean the lenscase.
- Last revision date: 2009/11

PRECAUTIONS:

- Once opened, the bottle should not be used after 90 days.
- Maintain the bottle well closed after use.
- Do not use the product if the safety seal is missing or broken.
- You may store your lenses up to 30 days in the closed lenscase
- Do not ingest the solution and keep away from children's reach.
- Keep at room temperature.
- Do not allow the dropper to come into contact with any surface to avoid contamination
- In case of irritation, consult your contact lens specialist.
- Before using any ophtalmic medicine, inform your doctor you are a contact lens user.
- Do not use if you are allergic to any of the components.
- Never use tap water to rinse your lenses or lens case.
- Do not use after the expiration date.
- Discard the solution in the lenscase after each use.

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