- EN MYLK CLEANER & CONDITIONER**
- Cleaning, disinfecting, rinsing, wetting, and storing for all types of rigid and gas permeable contact lenses.**
- INSTRUCCIONES DE USO:**
1. Wash your hands carefully with water and soap before handling your lenses.
 2. Place the lens on the palm of your hand.
 3. Add a few drops of MYLK Cleaner on the surface of each lens.
 4. Rub gently with the tip of your finger until a little foam appears.
 5. Rinse your lenses abundantly with MYLK Conditioner to eliminate all residue of the cleaning solution.
 6. Place your lenses in the lens case and cover them with MYLK Conditioner.
 7. To guarantee an adequate disinfection, lenses should be soaked in that solution for at least 4 hours or overnight.
 8. Before fitting your lenses, add a few drops of MYLK Conditioner on each lens.
 9. Keep your lens case clean and use a fresh supply of solution every day.

Last revision date: 2009/11

- PRECAUTIONS:**
- MYLK Cleaner causes irritation if not rinsed abundantly. The solution should never come in contact with your eyes.
 - Once opened, the bottles should not be used after 60 days.
 - Keep bottles well closed after use and keep away from children's reach.
 - You may store your lenses up to 30 days in the closed lens case.
 - Do not touch the dropper or dirty it.
 - In case of eye irritation consult a specialist.
 - Do not use these products with soft contact lenses.
 - Do not use the product if the safety seal is missing or broken.
 - Never use tap water to rinse your lenses or lens case.
 - Do not use after the expiration date.
 - Keep at room temperature.
 - Before using any ophthalmic medicine, inform your doctor you are a contact lens user.
 - Do not use if you are allergic to any of the components.
 - Discard the solution in the lens case after each use.

- COMPOSITION:**
- MYLK Cleaner: anionic detergent 10% in aqueous solution.
 - MYLK Conditioner: viscous agent, EDTA 0,1%, polyhexanide 0,0002% in an isotonic, buffered and sterile solution.

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