	CO.D. 227017	Rev. 01	Data 2024/07	LINEA Soleko
Prodotto PLI	Tipo LEAFLET			
FUSTELLA 200x450 mm	COLOR CMYK			

YALFRESH

MULTI-PURPOSE
CONTACT LENS SOLUTION

EN
PROPERTIES AND INDICATIONS: Yalfresh is a multipurpose solution containing hyaluronic acid, Clean, disinfects, lubricates, stores, wets, rinses and removes deposits. Suitable for soft contact lenses including coloured lenses. Yalfresh has an elevated antibacterial, cleaning and anti-protein action. The sodium hyaluronate gives the Yalfresh formula a unique balance, offering an effective antibacterial action, superior comfort and long-lasting hydration.

INSTRUCTIONS FOR USE: The vast majority of problems associated with wearing contact lenses are caused by failure to accurately follow the contact lens care instructions. To avoid problems and to ensure comfortable contact lens wear, it is important to follow steps 1, 2 and 3 every day.

Step 1 – Cleaning
Wash, rinse and dry your hands before handling the lenses. Place one lens in the palm of your hand and apply a few drops of Yalfresh onto each. Gently rub both sides of the lens with the tip of your forefinger for 15–20 seconds.

Step 2 – Rinsing
After cleaning, rinse your lenses with Yalfresh.

Step 3 – Disinfection
Remove the cap from your contact lens case. Place your left lens in the left compartment, and your right lens in the right compartment. Fill the case almost to the brim with Yalfresh and securely tighten the cap. Leave the lenses to soak for at least 4 hours or overnight. **PRECAUTIONS:** It is recommended that the lenses be re-disinfected after prolonged periods of storage. In seven days or longer, Discard contact lenses 60 days after opening. Do not use if tamper-evident seal is damaged. Adhere to the instructions recommended by your contact lens practitioner. Store between 5 and 30°C. Do not use the product after the expiry date.

WARNINGS: In the event of persistent eye irritation, discontinue use immediately. Remove the lens from the eye and consult your contact lens practitioner. Keep all contact lens care products away from children. Do not remove the cap from the case when the contact lenses are stored. Do not allow the nozzle tip to come into contact with any surfaces. Do not rinse the contact lenses or lens case directly with tap water. Never swallow the Yalfresh solution. Consult your contact lens practitioner before using any eye medication while wearing your contact lenses. Your lens storage case should be cleaned regularly and changed frequently as recommended by your contact lens practitioner. Any serious incidents that occur in relation to the device should be reported to Soleko S.p.A. and the competent authority of the Member State in which the user is established.

CONTRAINDICATIONS: Do not use the product if you are allergic or sensitive to any of the ingredients.

CONTENTS: PHMB 0.00015%, dibasic sodium phosphate, monobasic sodium phosphate, sodium chloride, disodium EDTA, sodium hyaluronate, hydroxyethylcellulose, purified water.

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KEY SYMBOLS USED ON THE PRIMARY CONTAINER, SECONDARY CONTAINER, INSTRUCTIONS FOR USE			
LOT	Manufacturer	MD	Symbol of Medical Device
LOT	Lot Number	MD	Consult Instructions for use
RECYCLE	Recycling symbol	CE	CE-marking
Use-by date	Date of expiration	Recycle	Recycling
Use-by date	Date of expiration	Recycle	Recycling
Use-by date	Date of expiration	Recycle	Recycling

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